

ORDERING DETAILS

410-010 Winner EVO CM2 2 Channel Combo

410-012 Winner EVO CM4 4 Channel Combo

STANDARD ACCESSORIES

410-160 5 cm² Applicator

GS5CM GelShot Therapeutic Ultrasound Coupling Discs - 5 cm (96 shots/box)

GSA5RM1 GelShot Adaptor - Richmar 5 cm

A1717A Lead Cord Set 1.8m (6 ft.)
CM2 (2 x); CM4 (4 x)

203-016 US Hospital Grade Power Cord
2.4m (8 ft.)

410-250 Instruction Manual

1 oz. (30 ml) AloeSound Lotion

1 pk. MultiStim Electrodes

OPTIONAL ACCESSORIES

410-040 StoreMore Cart

410-200 AutoSound Hands-Free Ultrasound
Applicator

410-115 Therapy Hammer 2 cm² / 5 cm² (2 in 1)

410-176 10 cm² Applicator (1 MHz)

600-123 Light Cluster Probe, Hands-Free
(500 mW - 870 nm LED / 633 nm SLD)

500-031 LaserPrism Laser Probe, Tethered
(200 mW / 785 nm)

500-030 LaserPrism Laser Probe, Cordless
(200 mW / 785 nm)

410-060A Laser Module - Winner EVO
(Required for Tethered Laser / Light
Probes)

K36031H5183BB-D UK Power Cord

K5A031H5183BB Euro Power Cord

TECHNICAL SPECIFICATIONS


Dimensions*: 12 1/2" x 16" x 9 1/2"
(31.75 x 40.64 x 24.13 cm)


Weight: CM2 - 22.5 lbs (10 kg)
CM4 - 24.5 lbs (11 kg)

Power: 115 / 230 VAC

Frequency: 60 / 50 Hz

Power Consumption: 110 Watts

Waveforms:  Interferential (4 Pole)

 Premodulated (2 Pole)

 High Volt

 Russian

 Microcurrent

LaserPrism Emitter - see page 9

*Dimensions are listed as L x W x H

WINNER EVO COMBO



The Winner EVO combination unit provides the clinician with the most comprehensive therapies in one unit with a reliable platform that is simple to use. Winner EVO CM4 can run 4 channels of electrical stimulation, ultrasound and laser &/or light all at the same time ensuring maximum utilization and return on investment in a busy clinic.

- Dual Frequency (1 & 3 MHz) Clinical Ultrasound with Industry Leading BNR
- Continuous or Pulsed Modes
- 5 cm², 2 / 5 cm² Therapy Hammer, 10 cm², or AutoSound Transducers Available
- Combination, Ultrasound, Electrotherapy, and/or Laser in One
- Tru-Stim Electrotherapy with Five Waveforms on-Board
- Therapy Profiles for Commonly Used Treatments
- Lead Cord Continuity Test on Board All Winner EVO Series
- Accessory Holders for Easy Applicator Organization
- Available AutoSound and Light Cluster Probe Allow for Hands-Free Therapy
- StoreMore Cart Is Available
- 3-Year Warranty
- Proudly Made in USA

 Visible and Invisible laser radiation. Avoid direct exposure to beam. Class IIIb laser product.