

ORDERING DETAILS

410-020	Therasound EVO (120V)
410-028	Therasound EVO (220V)
410-026	Therasound 3 Standard (not shown)

STANDARD ACCESSORIES

THERASOUND EVO

410-115	Therapy Hammer Ultrasound Applicator (2 & 5 cm ² in 1)
410-125	Instruction Manual
GS5CM	GelShot Therapeutic Ultrasound Coupling Discs - 5 cm
GSA5RM1	GelShot Adaptor - 5 cm
GSA2RM1	GelShot Adaptor - 2 cm
	AloeSound Lotion 1 oz. (30 ml)

THERASOUND 3 STANDARD

410-160	5 cm ² Applicator
410-125	Instruction Manual
GS5CM	GelShot Therapeutic Ultrasound Coupling Discs - 5 cm
GSA5RM1	GelShot Adaptor - Richmar 5 cm
	AloeSound Lotion 1 oz. (30 ml)

OPTIONAL ACCESSORIES

410-040	StoreMore Cart
410-042	Therasound Cradle (Mount unit to cart)
410-200	AutoSound Hands-Free Ultrasound Applicator
410-160	5cm ² Ultrasound Applicator
410-176	10cm ² Ultrasound Applicator (1MHz)
600-123	Light Cluster Probe, Hands-Free (500mW - 870nm LED/633nm SLD)
500-031	LaserPrism Laser Probe, Tethered (200 mW / 785 nm)
500-030	LaserPrism Laser Probe, Cordless (200 mW / 785 nm)
410-050A	Laser Module - Therasound EVO (Required for Tethered Laser / Light Probes)

TECHNICAL SPECIFICATIONS

Dimensions*:	9 1/2" x 12" x 6 1/4" (24.13 x 30.48 x 15.87 cm)
Weight:	5.6 lbs (2.5 kg)
Power:	115 VAC at 1 Amp / 230 VAC at .5 Amp
Frequency:	60 / 50 Hz

LaserPrism Probes - see page 8
*Dimensions are listed as L x W x H

US Patent No. 5558623


THERASOUND EVO



Richmar's heritage is ultrasound. The legacy products from Richmar were among the first ultrasound devices engineered more than forty years ago. That same knowledge, craftsmanship and quality reside in today's Therasound. A dual frequency clinical ultrasound device that has an industry leading BNR, Therasound can also become a hybrid ultrasound/laser modality with two independent channels. This field upgrade can be done now or later.

- Dual Frequency Clinical Ultrasound (1 & 3 MHz) with Industry Leading BNR
- Continuous or Pulsed Duty Cycles
- Patented Therapy Hammer 2 & 5 cm² Combined in One Transducer
- Therasound EVO Upgradeable to Accept Laser
- Both AutoSound and Light Cluster Probe Allow for Hands-Free Therapy
- With the Inclusion of a Therasound Cradle, Integration with a StoreMore Cart Is Available
- 3 Year Warranty
- Proudly Made in USA

5cm² Applicator 3MHz frequency
Acoustical Measurement @ 5ml
BNR: 2.87 ERA: 4.97
(FDA test standards)



5cm² Applicator 1MHz frequency
Acoustical Measurement @ 5ml
BNR: 2.94 ERA: 4.57
(FDA test standards)

