

## ORDERING DETAILS

410-200 AutoSound Hands-Free Ultrasound Applicator

## STANDARD ACCESSORIES

A1220 Swivel Attachment w/  
Three Fabrifoam Velcro Straps

400-046 AutoGel Pads (50x / case)

## TECHNICAL SPECIFICATIONS

Dimensions\*: 3.0" x 4.2" x 1.5"  
(7.62 x 10.7 x 3.8 cm)

Weight: 0.8 lbs (0.36 kg)

Power: 1 MHz = 0.0 - 7.5 W  
3 MHz = 0.0 - 5.0W

*\*Dimensions are listed as L x W x H*

## AUTOSOUND



Maximizing clinical efficiency while ensuring positive clinical outcomes can be achieved through the utilization of the time-tested AutoSound.

- Patented, Hands-Free Ultrasound Applicator
- 1 MHz, 3 MHz As Well As 1 & 3 MHz Sweep
- Applicator Includes Swivel Attachment with Three Fabrifoam Velcro Straps for Easy Fixation to the Body
- 4 Individual 3.5 cm<sup>2</sup> Transducers in Each AutoSound Applicator
- Sequences Between Each of the 4 3.5 cm<sup>2</sup> Transducers for 1 Second Each
- Easily Accessorize the Winner EVO CM2 / CM4 or Therasound EVO Device for Hands-Free Therapy
- 1 Year Warranty
- Proudly Made in USA

