

COMPREHENSIVE REHABILITATION & RESTORATIVE TRAINING

APT PERFORMER



APT PERFORMER Features & Benefits



Richmar's APT Performer is a multifaceted rehab training device for upper and lower extremities with simple-to-use features that benefit both patients and clinicians. It is used to improve functional outcomes by enhancing strength and coordination, ultimately restoring patients' ability to participate in their regular daily activities. The APT Performer includes industry-standard strength and symmetry training modes for upper and lower extremities, and is also [the first-and-only active-passive trainer on the market to incorporate fine motor skills training.](#)

ORDER DETAILS

DQAAPT APT Performer

TECHNICAL SPECIFICATIONS

Dimensions 28" x 30" x 42"
Rated Voltage AC 100-240 V
Frequency 50-60 Hz

*Fine motor training is available for upper extremities only

** Compatible with Polar H10 Sensors; not included

1297.20.09.B ©2020 Compass Health Brands Corp.



Symmetry Training

guides bilateral balance of muscle strength and coordination



Strength Training

aims to improve muscle conditioning with set force production goals



Endurance Training

encourages patients to increase their longevity of exercise output



ROM Training

maintains muscle and soft-tissue pliability, improves joint mobility, and helps prevent muscle atrophy



Fine Motor Training

measures hand-eye coordination on the basis of reaction time and accuracy



Biofeedback

prompts including auditory and visual feedback, on-screen treatment summary and HIPAA-compliant summary download



Intelligent Spasm Detection

with adjustable settings for spasm-recognition and remediation based on the Modified Ashworth Scale (MAS)



Compass Health Brands Corp.
Toll Free 1.888.549.4945
6753 Engle Road
Middleburg Heights, OH 44130
richmarweb.com