

PROFESSIONAL PAIN MANAGEMENT

THERATOUGH® CX4, THERATOUGH® CX2 & THERATOUGH® EX4





ELECTROTHERAPY ELECTROSTIMULATION

TheraTouch® EX4

- 8 clinical waveforms
 - Interferential
 - Premodulated
 - TENS (Symmetrical Biphasic & Asymmetrical Biphasic)
 - NMS
 - Russian
 - High Volt
 - Microcurrent
 - RAAS: New, proprietary protocol for muscle movement assistance
- Selectable evidence-based clinical protocols
- 4 channels of stimulation, use independently or in conjunction
- English and Spanish language options
- Save up to 99 user-defined favorites
- Patient interrupt switch
- Full-color touchscreen interface
- Optional therapy cart



COMBINATION THERAPY ELECTROSTIMULATION + ULTRASOUND

TheraTouch® CX4

- 5 cm² ultrasound applicator
 - 1 & 3 MHz ultrasound frequencies
- 8 clinical waveforms
 - Interferential
 - Premodulated
 - TENS (Symmetrical Biphasic & Asymmetrical Biphasic)
 - NMS
 - Russian
 - High Volt
 - Microcurrent
 - RAAS: New, Proprietary protocol for muscle movement assistance
- Selectable evidence-based clinical protocols
- 4 channels of stimulation, use independently or in conjunction
- English and Spanish language options
- Save up to 99 user-defined favorites
- Patient interrupt switch
- Full-color touchscreen interface
- Optional therapy cart

TheraTouch® CX2

NOW AVAILABLE — TheraTouch CX2 provides same features of TheraTouch CX4 in a more compact two-channel configuration.

FREE VIDEO TRAINING AVAILABLE

Learn the navigation and operation of these devices with our online video training series, including application suggestions using the pre-programmed clinical protocols. Preview it today, at: richmarweb.com/how-to-cx4



INDUSTRY'S BEST
3 YEAR
WARRANTY