

ORDERING DETAILS

A 600-123 Light Cluster Probe – Hands Free

B 500-029 LaserPrism Pro-1100 Console

410-060A Laser Module - Winner EVO

410-050A Laser Module - Therasound EVO

STANDARD ACCESSORIES

LASERPRISM PRO-1100 CONSOLE

MCT 503 Light Cluster Probe – Hands Free (2x)

LCT 201 Laserprism – Tethered

600-207 Laser goggles (2x)

TECHNICAL SPECIFICATIONS

LIGHT CLUSTER PROBE

Optical Output Power: 500 mW +/- 20%

Energy Delivered: 1.0 Joules / cm² in 45 sec

Wavelength: 870 nm

Irradiance: 22.2 mW / cm² +/-20%

Duty Cycle: 100% continuous mode

Treatment area: 22.48 cm²

Laser Safety Classification: Class IIIB

LASERPRISM PRO-1100 CONSOLE

Dimensions: 19" x 17" x 5.5"
(48 x 43 x 14 cm)

Weight: 5.5 lbs. (2.5 kg)

Power: 100 / 240 VAC,

Frequency: 50 / 60Hz

Duty Cycle: 100% continuous mode
or 80% for 3 pulse mode
switch mode power
supply

**Dimensions are listed as L x W x H*

Note: Laser module required for light cluster probe connection to a Winner EVO or Therasound EVO



A



B

Light Cluster Probe

- Hands-Free SLD Light Therapy
- 500 mW Total Power, 870 nm LED / 633 nm SLD Frequency
- Delivers 1 Joule / cm² Every 45 Seconds
- Treatment Time Can Be Adjusted If Used in Conjunction with the Tethered Laser
- 1 Year Warranty

LaserPrism Pro-1100 Console

- Equipped with 2x Hands-Free Light Cluster Probes and 1x LaserPrism Emitter
- Includes 785 nm Infrared GaAlAs Laser Diode Applicator with 200 mW Total Power
- Includes LED Visible Red (633 nm) and SLD Infrared (870 nm) Hands Free Light Cluster with 500 mW of Total Power
- Continuous or Pulsed Wave Modes
- 3 Port Console for LLLT, LED/SLD Treatments
- Power Console with Built-in Treatment Guidelines Provide Parameter Selection
- Automatically Calculates and Displays Treatment Times
- 1 Year Warranty

Laser Module (not pictured)

- Optional Upgrade to Any Winner EVO Allows Dual Operation of Light/Laser While Performing Up to 4 Channels of Electrotherapy All Simultaneously
- Optional Upgrade to Therasound EVO Allows Dual Operation of Light/Laser and Ultrasound Simultaneously
- Two Ports for LaserPrism and/or Light Cluster Plug-in
- 1 Year Warranty