

## ORDERING DETAILS

HT-R12-S	<b>HydraTherm Standard 120V</b>
HT-R12-S220	<b>HydraTherm Standard 220V</b> Without Packs, Divider System
HT-R12-SW	<b>HydraTherm Standard 120V</b>
HT-R12-SW220	<b>HydraTherm Standard 220V</b> Includes 3 Standard, 3 Universal, 3 Cervical, & 3 Oversize HydraHeat Packs, Divider System
HT-R12-SR	<b>HydraTherm Standard 120V</b>
HT-R12-SR220	<b>HydraTherm Standard 220V</b> Without Packs, Rack System
HT-R12-SRW	<b>HydraTherm Standard 120V</b>
HT-R12-SRW220	<b>HydraTherm Standard 220V</b> Includes 3 Standard, 3 Universal, 3 Cervical, & 3 Oversize HydraHeat Packs, Rack System
HT-R12-D	<b>HydraTherm Deluxe 120V</b>
HT-R12-D220	<b>HydraTherm Deluxe 220V</b> Includes Drainage Pump, Without Packs; Divider System
HT-R12-DW	<b>HydraTherm Deluxe 120V</b>
HT-R12-DW220	<b>HydraTherm Deluxe 220V</b> Includes Drainage Pump, Includes 3 Standard, 3 Universal, 3 Cervical, & 3 Oversize HydraHeat Packs, Divider System
HT-R12-DR	<b>HydraTherm Deluxe 120V</b>
HT-R12-DR220	<b>HydraTherm Deluxe 220V</b> Includes Drainage Pump, Without Packs, Rack System
HT-R12-DRW	<b>HydraTherm Deluxe 120V</b>
HT-R12-DRW220	<b>HydraTherm Deluxe 220V</b> Includes Drainage Pump, Includes 3 Standard, 3 Universal, 3 Cervical, & 3 Oversize HydraHeat Packs, Rack System

## OPTIONAL ACCESSORIES

HT-R12-RCK	Rack System
HT-R12-DIV	Divider System (includes 3 clip-on shelves)
HTR121	Fold-Away Side Shelf
HT-R12-DIVS	Clip-On Shelves

## TECHNICAL SPECIFICATIONS

Dimensions:	30" W x 20" D x 33" H (76 x 50 x 84 cm)
Weight:	Standard - 70 lbs (31 kg) Deluxe - 75 lbs (34 kg)
Power:	110 - 120 VAC / 220 - 240 VAC
Frequency:	50 / 60 Hz
Power Consumption:	1000 Watts
Adjustable Temp. Range:	120 °F - 160 °F (71 °C - 49 °C)
Safety thermal cut-out Temp.:	175 °F +/- 5 °F (79 °C +/- 3 °C)
Digital Temp. Control Accuracy:	+/- 5%
Water capacity:	18 Gal (68 L) (without heat packs)

## HYDRATHERM



HydraTherm with Divider System

HydraTherm heating units set the standard for modern moist heat therapy with an updated modern design, featuring a rust-free composite tank that will not dent, and unique hot pack storage with a 12 pack capacity.

- Only 2 Recommended Cleanings Per Year Thanks to Rust-Free Composite Tank, Compared to 24 for Traditional Stainless Tanks
- Digital Thermostat with Easy to Read Temperature Display Gives Therapist Ability to Precisely Target Treatment Temperature, Increasing Patient Safety
- Safe Operation with Low Water Cut-Off, and 24 Volt Control Unit, Compared to 110 Volt Competitors
- HydraTherm Deluxe Models Include an Electric Drainage Pump and Fold-Away Side Shelf
- Hot Pack Storage Options Include the Divider System That Provides Tong-Free Access to Hot Packs, or the Rack System That Allows for Additional Pack Storage on Top of Racks.
- Use in Conjunction with HydraHeat Packs for a Clean, Hygienic Moist Heat Therapy System
- 1 Year Warranty
- Proudly Made in USA



HydraTherm with Rack System