



ORDERING DETAILS

HP-1713-OS	Oversize 17" x 13" (43 x 33cm)
HP-1905-UN	Universal 19" x 5" (48 x 13cm)
HP-1711-CX	Neck Contour/Cervical 17" x 11" (43 x 28cm)
HP-1711-ST	Standard 17" x 11" (43 x 28cm)

HydraHeat Packs paired with HydraTherm heating units create a complete moist heat therapy system that provides patient comfort, cleanliness and a time-efficient clinic.

- The Polymer Skin of HydraHeat Packs Are Easily Disinfected with Antibacterial Wipes, Helping Reduce Cross Contamination Between Patients
- HydraHeat Packs Do Not Require High Temperatures to Inhibit Bacterial Growth As with Clay Packs, So Lower, Safer Temperatures Within the Therapeutic Range Can Be Used to Treat Patients, Minimizing Burn Risk
- Lower Treatment Temperatures Allow HydraHeat Packs to Be Covered with Only a Towel or Pillow Case, Not Bulky, Expensive, and Unsanitary Terry Covers
- The Sureseam Process Used to Seal Each Pack Lessens the Chances of Bursting or Leaking its Contents Into the Tank, or Onto the Patient
- Polymer Skin in Combination with the Non-Toxic Gel Interior Make for a Pliable Pack That Conforms to Various Treatment Areas and Allows a Greater Area of Contact for Heat Transfer
- HydraHeat Packs Maintain Heat Within the Therapeutic Range for Up to 40 Minutes, and Reheat to Treatment Temperature Within 15 Minutes
- Not Made with Natural Rubber Latex