

superstim
PREMIUM ELECTRODES

— & —

multistim
ELECTRODES



Another innovative idea from


 **Richmar**[®]
Engineering Tomorrow's
Healthcare

superstim

PREMIUM ELECTRODES



superstim premium
Silver carbon
Electrodes



Richmar electrode technology

Over 40 years of electrotherapy manufacturing experience, we have developed a line of superior quality electrodes. Made in the U.S.A. and made with a proprietary bonding and moisture barrier technology that eliminates worries when applying moist heat therapy. It also provides superior performance in the industry. The electrodes are constructed with two hydrogel layers that provide comfortable and uniform stimulation.

multistim

ELECTRODES



multistim carbon
Film electrodes

superstim & multistim by Richmar

- Available in cloth or foam
- Strongest lead wire pull strength in the industry
- Premium dual layer 35 mil gel
- Gel formulation produces maximum dry-out resistance
- Optimum distribution of electrical current
- Proudly made in the U.S.A.

Richmar's experience has provided Richmar the insight to create SuperStim. In the U.S.A., the SuperStim and MultiStim Electrodes are made with a top layer called Leadstay. This eliminates delamination and provides the highest lead wire tensile strength in the industry. The electrodes feature drogel layers of premium 35 mil thickness and

Ordering Details

REF	Description
superstim Silver-Carbon Electrodes	
400-895	5cm (2") Round (Cloth)
400-894	5cm (2") Round (Foam)
400-899	5cm x 5cm (2" x 2") Square (Cloth)
400-898	5cm x 5cm (2" x 2") Square (Foam)
400-893	7.6cm (3") Round (Cloth)
400-892	7.6cm (3") Round (Foam)
400-864	3cm (1.25") Round (Cloth)
400-869	3cm (1.25") Round (Foam)
400-858	4cm x 9cm (1.5" x 3.5") Rectangle (Cloth)
400-857	4cm x 9cm (1.5" x 3.5") Rectangle (Foam)
400-889	5cm x 9cm (2" x 3.5") Rectangle (Cloth)
400-888	5cm x 9cm (2" x 3.5") Rectangle (Foam)
400-854	7.6cm x 13cm (3" x 5") Rectangle (Cloth)
400-853	7.6cm x 13cm (3" x 5") Rectangle (Foam)

REF	Description
multistim Carbon Electrodes	
400-872	5cm (2") Round (Cloth)
400-873	5cm (2") Round (Foam)
400-877	5cm x 5cm (2" x 2") Square (Cloth)
400-878	5cm x 5cm (2" x 2") Square (Foam)
400-881	7.6cm (3") Round (Cloth)
400-882	7.6cm (3") Round (Foam)
400-863	3cm (1.25") Round (Cloth)
400-862	3cm (1.25") Round (Foam)
400-856	4cm x 9cm (1.5" x 3.5") Rectangle (Cloth)
400-855	4cm x 9cm (1.5" x 3.5") Rectangle (Foam)
400-879	5cm x 9cm (2" x 3.5") Rectangle (Cloth)